

Streptavidin Monoclonal Antibody For Flow Cytometry - Sa V1, Sa V2

Item Number: CM0012001

Introduction

Streptavidin is a monoclonal antibody with species-independent cross-reactivity, ideal for flow cytometry applications. Sourced as an IVD raw material for diagnostic developers and researchers.

[Learn More](#)

Parameter	Value
Product Name	Streptavidin
Alias	SA V1, SA V2
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Species independent
Shipping	Ice bag