

Stx1A Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - Q16623

Item Number: CM0002037

Introduction

High-quality STX1A rabbit monoclonal antibody validated for WB, IF-P, IHC-P, and ELISA. Cross-reacts with human, mouse, rat, targeting the neuronal SNARE protein syntaxin-1A (UniProt Q16623).

[Learn More](#)

Parameter	Value
Product Name	STX1A Rabbit mAb
Remarks/Alias	STX1; HPC-1; P35-1; SYN1A; STX1A
Species	Human
GeneID (Human)	6804
GeneID	6804
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human STX1A (Q16623)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16623
Protein Weight	33kDa
Shipping	Ice bag