

Sufu Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q9Umx1

Item Number: CM0013348



Introduction

SUFU rabbit polyclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Recognizes human, mouse, and rat SUFU, a negative regulator of the Hedgehog pathway. Ideal for cancer and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	SUFU Rabbit pAb
Remarks/Alias	SUFUH; JBTS32; SUFUXL; PRO1280; SUFU
Species	Human
GeneID (Human)	51684
GeneID	51684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 325-484 of human SUFU (NP_057253.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UMX1
Protein Weight	54kDa
Shipping	Ice bag