


Sumo2/3 Rabbit Polyclonal Antibody (Wb, Elisa) - P61956 / P55854

Item Number: CM0022476



Introduction

Rabbit polyclonal antibody against human SUMO2/3, validated for Western blot and ELISA. Detects SUMO2 and SUMO3 cross-reactive with mouse and rat; suitable for studying sumoylation and proteasomal degradation.

[Learn More](#)

Parameter	Value
Product Name	SUMO2/3 Rabbit pAb
Remarks/Alias	HSMT3; SMT3B; SUMO3; Smt3A; SMT3H2; SUMO2/3
Species	Human
GeneID (Human)	6613/6612
GeneID	6613/6612
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human SUMO2/3 (NP_008868.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P61956/P55854
Protein Weight	11 kDa
Shipping	Ice bag