

Surf4 Rabbit Pab - O15260

Item Number: CM0020196

Introduction

SURF4 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, IF/ICC, ELISA. It detects human SURF4 (also known as ERV29) and cross-reacts with mouse and rat. SURF4 functions as an endoplasmic reticulum cargo receptor for lipoprotein secretion.

[Learn More](#)

Parameter	Value
Product Name	SURF4 Rabbit pAb
Remarks/Alias	ERV29; SURF4
Species	Human
GeneID (Human)	6836
GeneID	6836
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15260
Protein Weight	30kDa
Shipping	Ice bag