

Symmetric Dimethyl-Histone H3-R17 Rabbit Monoclonal Antibody (Q16695/P68431) For Wb, Db, Ihc-P, Elisa

Item Number: CM0000332



Introduction

Rabbit monoclonal antibody targeting symmetric dimethylation at R17 of human Histone H3. Validated for Western Blot, Dot Blot, IHC-P, and ELISA. Cross-reacts with human, mouse, rat, and other species. Ideal for epigenetic research and histone modification analysis.

[Learn More](#)

Parameter	Value
Product Name	Symmetric DiMethyl-Histone H3-R17 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Symmetric DiMethyl-Histone H3-R17
Species	Human
GeneID (Human)	8290/ 8350
GeneID	8290/8350
Immunogen	A synthetic symmetric dimethylated peptide around R17 of human Histone H3-R17 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag