

Symmetric Dimethyl-Histone H3-R2 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Ihc, If, Elisa, Chip-Seq

Item Number: CM0000342



Introduction

Rabbit polyclonal antibody against symmetric dimethylated arginine 2 of histone H3 (Q16695/P68431). Validated for WB, IHC-P, IF/ICC, ELISA, and ChIP-seq; reacts with human, mouse, rat, and predicted wide range. Ideal for epigenetic research and histone modification studies.

[Learn More](#)

Parameter	Value
Product Name	Symmetric DiMethyl-Histone H3-R2 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Symmetric DiMethyl-Histone H3-R2
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic symmetric dimethylated peptide around R2 of human histone H3 (NP_003520.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag