

Symmetric Dimethyl-Histone H3-R26 Rabbit Pab (Q16695/P68431) For Wb, If/Icc, Elisa

Item Number: CM0000290



Introduction

Rabbit polyclonal antibody detecting Symmetric DiMethyl-Histone H3-R26, suitable for WB, IF/ICC, ELISA applications. Cross-reactivity with human, mouse, and rat. Ideal for epigenetic research and IVD development.

[Learn More](#)

Parameter	Specification
Product Name	Symmetric DiMethyl-Histone H3-R26 Rabbit pAb
Remarks/Alias	H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3Fj; H3C10; H3C11; HIST1H3j; Symmetric DiMethyl-Histone H3-R26
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	A synthetic symmetric dimethylated peptide around R26 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag