

Symmetric Dimethyl-Histone H3-R8 (Q16695/P68431) Rabbit Mab For Western Blot, Ihc, And If/Icc

Item Number: CM0000288



Introduction

Symmetric DiMethyl-Histone H3-R8 Rabbit mAb validated for WB, DB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat and wide-range species. Detects a 16kDa band; suitable for histone modification and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Symmetric DiMethyl-Histone H3-R8 Rabbit mAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Symmetric DiMethyl-Histone H3-R8
Species	Human
GeneID (human)	8290/8350
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic symmetric dimethylated peptide around R8 of human Histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	16kDa
Shipping	Ice bag