


Symmetric Dimethyl-Histone H3-R8 Rabbit Pab (Q16695/P68431) For Wb, Ihc-P, If/Icc, Elisa, Chip-Seq

Item Number: CM0000372



Introduction

Rabbit polyclonal antibody specific for symmetric dimethylation at R8 of human histone H3 (UniProt Q16695/P68431). Validated for WB, IHC-P, IF/ICC, ELISA, and ChIP-seq. Cross-reacts with human, mouse, rat, and other species.

[Learn More](#)

Parameter	Specification
Product Name	Symmetric DiMethyl-Histone H3-R8 Rabbit pAb
Remarks/Alias	H3t; H3.4; H3/g; H3FT; H3C16; HIST3H3; Symmetric DiMethyl-Histone H3-R8
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide - A synthetic symmetric dimethylated peptide around R8 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA, ChIP-seq
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15 kDa
Shipping	Ice bag