

Symmetric Dimethyl-Histone H3-R8 Rabbit Polyclonal Antibody (Q16695/P68431) For Wb, Elisa, Chip

Item Number: CM0000349



Introduction

Rabbit polyclonal antibody against symmetric dimethyl-Histone H3-R8 (SWISS Q16695/P68431). Suitable for WB, ELISA, and ChIP. Recognizes human, mouse, rat, and other species (wide range predicted). Ideal for epigenetic and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	Symmetric DiMethyl-Histone H3-R8 Rabbit pAb
Remarks/Alias	H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Symmetric DiMethyl-Histone H3-R8
Species	Human
GeneID (human)	H3
GeneID	8290/8350
Immunogen	Synthetic peptide: A synthetic symmetric dimethylated peptide around R8 of human histone H3 (NP_003520.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	Q16695/P68431
Protein Weight	15kDa
Shipping	Ice bag