

Symmetric Dimethyl-Histone H4-R3 (P62805) Rabbit Polyclonal Antibody For Wb, Ihc, If/Icc, Elisa

Item Number: CM0000302



Introduction

Symmetric DiMethyl-Histone H4-R3 Rabbit pAb (P62805) validated for WB, DB, IHC-P, IF/ICC, ELISA. Recognizes symmetric dimethylation at Arg3 of Histone H4. Reacts with human, mouse, rat and predicted wide range. Ideal for epigenetic research and histone modification analysis.

[Learn More](#)

Parameter	Value
Product Name	Symmetric DiMethyl-Histone H4-R3 Rabbit pAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A; Symmetric DiMethyl-Histone H4-R3
Species	Human
Gene Symbol (Human)	H4
Gene ID	8359
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Histone H4 (NP_003529.1).
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag