


Symmetric Dimethyl-Histone H4-R3 Rabbit Monoclonal Antibody (P62805) For Epigenetics Research

Item Number: CM0000251



Introduction

Rabbit monoclonal antibody specific for symmetric dimethyl-Histone H4 at Arg3 (P62805). Validated for Western blotting, IF/ICC, dot blot, and ELISA. Cross-reacts with human and rat histone H4. Ideal for studying histone code and epigenetic regulation mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Symmetric DiMethyl-Histone H4-R3 Rabbit mAb
Remarks/Alias	H4; H4/n; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4F2; H4FN; FO108; H4-16; H4C11; H4C12; H4C13; H4C15; H4C16; HIST2H4; HIST2H4A
Species	Human
GeneID (human)	8359
GeneID	8359
Immunogen	Synthetic peptide corresponding to amino acids 1-100 of human Histone H4 (NP_003539.1)
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IF/ICC, DB, ELISA
Cross Reactivity	Human, Rat
SWISS	P62805
Protein Weight	11kDa
Shipping	Ice bag