

Syntaxin 4 Rabbit Mab - Q12846

Item Number: CM0006231

Introduction

High-quality Syntaxin 4 Rabbit monoclonal antibody for WB and ELISA applications. Cross-reacts with human, mouse, and rat samples. Ideal for studying SNARE-mediated vesicle trafficking and glucose transport.

[Learn More](#)

Parameter	Value
Product Name	Syntaxin 4 Rabbit mAb
Remarks/Alias	STX4A; p35-2; Syntaxin 4
Species	Human
GeneID (Human)	6810
GeneID	6810
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human Syntaxin 4 (NP_004595.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12846
Protein Weight	34 kDa
Shipping	Ice bag