


Syvn1 Rabbit Polyclonal Antibody For Wb, Elisa - Q86Tm6

Item Number: CM0025244



Introduction

Rabbit polyclonal antibody targeting human SYVN1 (HRD1), an E3 ubiquitin ligase involved in ERAD. Validated for WB and ELISA; cross-reacts with mouse. Useful for studying ER stress, apoptosis, and ubiquitination pathways.

[Learn More](#)

Parameter	Value
Product Name	SYVN1 Rabbit pAb
Remarks/Alias	DER3; HRD1; SYVN1
Species	Human
GeneID (Human)	84447
GeneID	84447
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 500-600 of human SYVN1 (NP_757385.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q86TM6
Protein Weight	68kDa
Shipping	Ice bag