

## T-Bet/Tbx21 Rabbit Mab For Wb, Ihc-P, Fc (Intra), Elisa - Q9UI17

Item Number: CM0004604

### Introduction

Recombinant rabbit monoclonal antibody targeting human T-bet/Tbx21 (TBX21; T-PET; T-bet) for WB, IHC-P, FC (intra), ELISA. Ideal for Th1 cell lineage and transcription factor research.

[Learn More](#)

| Parameter        | Value                                                                                                                   |
|------------------|-------------------------------------------------------------------------------------------------------------------------|
| Product Name     | T-bet/Tbx21 Rabbit mAb                                                                                                  |
| Remarks/Alias    | TBX21; T-PET; T-bet; TBET; TBLYM; T-box 21                                                                              |
| Species          | Human                                                                                                                   |
| GeneID (Human)   | 30009                                                                                                                   |
| GeneID           | 30009                                                                                                                   |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 322-535 of human T-bet/Tbx21(NP_037483.1) |
| Source           | Rabbit                                                                                                                  |
| Category         | Monoclonal Antibodies                                                                                                   |
| Application      | WB, IHC-P, FC (intra), ELISA                                                                                            |
| Cross Reactivity | Human                                                                                                                   |
| SWISS            | Q9UL17                                                                                                                  |
| Protein Weight   | 58kDa                                                                                                                   |
| Shipping         | Ice bag                                                                                                                 |