

Taf15 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q92804

Item Number: CM0011523



Introduction

TAF15 rabbit monoclonal antibody validated for Western blot, IHC-P, IF/ICC, IP, and ELISA. Cross-reacts with human, mouse, and rat. Ideal for transcription and RNA-binding protein studies. SWISS: Q92804.

[Learn More](#)

Parameter	Value
Product Name	TAF15 Rabbit mAb
Remarks/Alias	Npl3, RBP56, TAF2N, TAFII68, TAF15
Species	Human
GeneID	8148
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92804
Protein Weight	62kDa
Shipping	Ice bag