


Taf6L Rabbit Polyclonal Antibody For Wb, Ip, Elisa - Q9Y6J9

Item Number: CM0016048



Introduction

Rabbit polyclonal antibody against TAF6L (PAF65A), validated for WB, IP, ELISA in human, mouse, rat. Suitable for studying PCAF complex, histone acetylation, and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	TAF6L Rabbit pAb
Remarks/Alias	PAF65A; TAF6L
Species	Human
Gene ID	10629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human TAF6L (NP_006464.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y6J9
Protein Weight	68kDa
Shipping	Ice bag