


Taf9 Rabbit Polyclonal Antibody For Wb And Elisa - Q16594

Item Number: CM0027665



Introduction

Rabbit polyclonal antibody against human TAF9 (TFIID subunit). Reacts with human and mouse. Suitable for WB and ELISA. Ideal for studies on transcription initiation and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	TAF9 Rabbit pAb
Remarks/Alias	TAF2G; TAFII31; TAFII32; MGC:5067; TAFII-31; TAFII-32; TAFIID32; STAF31/32
Species	Human
GeneID (human)	6880
GeneID	6880
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q16594
Protein Weight	29kDa
Shipping	Ice bag