

Tarbp2 Rabbit Pab For Wb - Q15633

Item Number: CM0027388



Introduction

TARBP2 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human and mouse TARBP2, a RISC-loading complex subunit critical for miRNA processing and RNA silencing. Suitable for viral infection studies.

[Learn More](#)

Parameter	Specification
Product Name	TARBP2 Rabbit pAb
Remarks/Alias	LOQS; TRBP; TRBP1; TRBP2; TARBP2
Species	Human
GeneID (human)	6895
GeneID	6895
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-366 of human TARBP2 (NP_599150.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15633
Protein Weight	39kDa
Shipping	Ice bag