


Taz Rabbit Monoclonal Antibody For Wb And Elisa - Q9Gzv5

Item Number: CM0006986



Introduction

TAZ rabbit monoclonal antibody validated for WB and ELISA, detects WWTR1/TAZ across human, mouse, and rat. Ideal for Hippo signaling, organ size control, and tumor suppression studies.

[Learn More](#)

Parameter	Value
Product Name	TAZ Rabbit mAb
Remarks/Alias	TAZ
Species	Human
GeneID (Human)	25937
GeneID	25937
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human TAZ (NP_056287.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9GZV5
Protein Weight	44kDa
Shipping	Ice bag