

Tbca Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 075347

Item Number: CM0015719



Introduction

TBCA Rabbit pAb validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat TBCA (tubulin-folding cofactor A), a 13 kDa chaperone essential for tubulin assembly and cytoskeleton integrity.

[Learn More](#)

Parameter	Specification
Product Name	TBCA Rabbit pAb
Remarks/Alias	TBCA
Species	Human
GeneID (human)	6902
GeneID	6902
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human TBCA (NP_004598.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075347
Protein Weight	13kDa
Shipping	Ice bag