

Tbx2 Rabbit Pab - Q13207

Item Number: CM0015263

Introduction

TBX2 rabbit polyclonal antibody suitable for Western blot and ELISA. Detects human, mouse, and rat TBX2, a transcription factor in development and EMT. Useful for cancer research.

[Learn More](#)

Parameter	Value
Product Name	TBX2 Rabbit pAb
Remarks/Alias	VETD; TBX2
Species	Human
GeneID (Human)	6909
GeneID	6909
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 604-712 of human TBX2 (NP_005985.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13207
Protein Weight	75kDa
Shipping	Ice bag