


Tbx5 Rabbit Polyclonal Antibody For Wb, Elisa - Q99593

Item Number: CM0014911



Introduction

Rabbit polyclonal antibody to TBX5, a T-box transcription factor crucial for heart development and limb formation. Suitable for WB and ELISA, it recognizes human, mouse, and rat TBX5. Ideal for developmental biology and cardiac research.

[Learn More](#)

Parameter	Value
Product Name	TBX5 Rabbit pAb
Remarks/Alias	HOS; TBX5
Species	Human
Gene ID (Human)	6910
Gene ID	6910
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human TBX5 (NP_852259.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99593
Protein Weight	58kDa
Shipping	Ice bag