

Tceb2 Rabbit Pab - Q15370

Item Number: CM0025131



Introduction

Rabbit polyclonal antibody against human TCEB2/Elongin-B. Suitable for WB, IP, ELISA. Cross-reacts with human and mouse. Ideal for transcription elongation and ubiquitin ligase research.

[Learn More](#)

Parameter	Value
Product Name	TCEB2 Rabbit pAb
Remarks/Alias	SIII; TCEB2
Species	Human
GeneID (Human)	6923
GeneID	6923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-118 of human TCEB2 (NP_009039.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15370
Protein Weight	13kDa
Shipping	Ice bag