

Tchh Rabbit Polyclonal Antibody For Wb And Elisa - Q07283

Item Number: CM0022787

Introduction

TCHH rabbit polyclonal antibody validated for WB and ELISA. Detects human and mouse Trichohyalin, a 254kDa intermediate filament-associated protein crucial for hair follicle inner root sheath and epidermal differentiation. Optimal for dermatology and keratinization studies.

[Learn More](#)

Parameter	Value
Product Name	TCHH Rabbit pAb
Remarks/Alias	THH; THL; TRHY; UHS3; TCHH
Species	Human
GeneID (Human)	7062
GeneID	7062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-300 of human TCHH (NP_009044.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q07283
Protein Weight	254kDa
Shipping	Ice bag