


Tet3 Rabbit Monoclonal Antibody For Wb And Elisa - Q8Bg87

Item Number: CM0004191



Introduction

CamelBio offers TET3 Rabbit Monoclonal Antibody validated for WB and ELISA in mouse. Studies epigenetic reprogramming and DNA demethylation in early embryonic development - Q8BG87.

[Learn More](#)

Parameter	Value
Product Name	TET3 Rabbit mAb
Remarks/Alias	B430006D22Rik; D230004J03Rik; TET3
Species	Mouse
GeneID (Human)	200424
GeneID	194388
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of mouse TET3 (Q8BG87)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8BG87
Protein Weight	195kDa
Shipping	Ice bag