


Tff1/Ps2 Rabbit Monoclonal Antibody For Wb, Ip, If-P, Elisa - P04155

Item Number: CM0006823



Introduction

Rabbit monoclonal antibody targeting human TFF1/pS2 (trefoil factor 1), validated for WB, IP, IF-P, ELISA. Cross-reacts with mouse and rat. Ideal for gastrointestinal and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	TFF1/pS2 Rabbit mAb
Remarks/Alias	pS2; BCEI; HPS2; HP1.A; pNR-2; D21S21; TFF1/pS2
Species	Human
GeneID (Human)	7031
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 25-84 of human TFF1/pS2. (NP_003216.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04155
Protein Weight	9kDa
Shipping	Ice bag