

Tfpi Rabbit Polyclonal Antibody For Wb, Elisa - P10646

Item Number: CM0023491



Introduction

High-quality TFPI Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, and rat. Targets tissue factor pathway inhibitor involved in blood coagulation. Ideal for research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	TFPI Rabbit pAb
Remarks/Alias	EPI, TFI, LACI, TFPI1, TFPI
Species	Human
GeneID (Human)	7035
GeneID	7035
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-224 of human TFPI (NP_001027452.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10646
Protein Weight	35kDa
Shipping	Ice bag