

Tfpt Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P0C1Z6

Item Number: CM0024885

Introduction

High-quality TFPT rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Targets human TFPT (INO80F), a component of the chromatin remodeling INO80 complex, promoting p53-independent apoptosis. Ideal for apoptosis and chromatin research.

[Learn More](#)

Parameter	Value
Product Name	TFPT Rabbit pAb
Remarks/Alias	FB1; amida; INO80F; TFPT
Species	Human
GeneID (Human)	29844
GeneID	29844
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-253 of human TFPT (NP_037474.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P0C1Z6
Protein Weight	28kDa
Shipping	Ice bag