

Tgf Beta 1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P01137

Item Number: CM0014547



Introduction

Rabbit polyclonal antibody against human TGF beta 1 (UniProt P01137). Validated for WB, IHC-P, IF/ICC, ELISA. Recognizes human, mouse, and rat TGF beta 1. Ideal for studying TGFβ signaling, fibrosis, and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	TGF beta 1 Rabbit pAb
Remarks/Alias	CED; LAP; DPD1; TGFB; IBDIMDE; TGFbeta; TGF-beta1; TGF beta 1
Species	Human
GeneID (Human)	7040
GeneID	7040
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-280 of human TGF beta 1 (NP_000651.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01137
Protein Weight	44 kDa
Shipping	Ice bag