


Tgf Beta Receptor Ii Rabbit Pab - P37173

Item Number: CM0017386



Introduction

Rabbit polyclonal antibody against TGF beta Receptor II (TGFB2) for Western blot, IF/ICC, and ELISA. Detects human, mouse, and rat targets. Suitable for TGF-beta signaling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	TGF beta Receptor II Rabbit pAb
Remarks/Alias	AAT3; FAA3; LDS2; MFS2; RIIC; LDS1B; LDS2B; TAAD2; TBR1; TBR-ii; TGFR-2; tbetaR-II; TGFbeta-RII; TGF beta Receptor II (TGFB2)
Species	Human
GeneID (Human)	7048
GeneID	7048
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-280 of human TGF beta Receptor II (TGFB2) (NP_003233.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P37173
Protein Weight	65 kDa
Shipping	Ice bag