

Tgfb2 Rabbit Pab For Wb, Elisa - P61812

Item Number: CM0027142



Introduction

Rabbit polyclonal antibody against human TGFB2, validated for WB and ELISA. Cross-reacts with mouse. Recombinant immunogen (aa 21-414). Targets Transforming growth factor beta-2 proprotein, a multifunctional cytokine regulating angiogenesis and heart development.

[Learn More](#)

Parameter	Value
Product Name	TGFB2 Rabbit pAb
Remarks/Alias	LDS4; G-TSF; TGF-beta2
Species	Human
GeneID (Human)	7042
GeneID	7042
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-414 of human TGFB2 (NP_003229.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	P61812
Protein Weight	48kDa
Shipping	Ice bag