


Thoc4/Allyref Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q86V81

Item Number: CM0009942



Introduction

THOC4/ALYREF Rabbit Monoclonal Antibody validated for Western blot, IHC-P, and ELISA. Detects human, mouse, and rat THO complex subunit 4, involved in mRNA export and processing. Ideal for cancer and RNA biology research.

[Learn More](#)

Specification	Details
Product Name	THOC4/ALYREF Rabbit mAb
Remarks/Alias	ALY; BEF; REF; THOC4; ALY/REF; THOC4/ALYREF
Species	Human
GeneID (Human)	10189
GeneID	10189
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 158-257 of human THOC4/ALYREF (Q86V81).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86V81
Protein Weight	27kDa
Shipping	Ice bag