

Thromboxane Synthase Rabbit Polyclonal Antibody For Wb, Elisa - P24557

Item Number: CM0015931

Introduction

Rabbit polyclonal antibody targeting thromboxane-A synthase (TBXAS1/CYP5A1) for WB and ELISA in human and mouse samples. Involved in thromboxane A2 synthesis and platelet aggregation. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Thromboxane synthase Rabbit pAb
Remarks/Alias	TS, TXS, CYP5, THAS, TXAS, CYP5A1, GHOSAL, BDPLT14, Thromboxane synthase
Species	Human
Gene ID (Human)	6916
Gene ID	6916
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 355-534 of human Thromboxane synthase (NP_001052.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P24557
Protein Weight	61kDa
Shipping	Ice bag