

Thyroid Peroxidase Rabbit Monoclonal Antibody For Ihc-P/Elisa - P07202

Item Number: CM0002702

Introduction

Thyroid Peroxidase Rabbit Monoclonal Antibody validated for IHC-P and ELISA. Detects human, mouse, and rat TPO. Ideal for thyroid hormone biosynthesis research and diagnostic assay development.

[Learn More](#)

Parameter	Description
Product Name	Thyroid Peroxidase Rabbit mAb
Remarks/Alias	MSA; TPX; TDH2A; Thyroid Peroxidase
Species	Human
Gene ID	7173
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Thyroid Peroxidase (P07202)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07202
Protein Weight	103kDa
Shipping	Ice bag