


Tiam1 Rabbit Pab For Wb, Ihc-P, Elisa - Q13009

Item Number: CM0014550



Introduction

TIAM1 Rabbit pAb is a polyclonal antibody validated for WB, IHC-P, and ELISA. It detects human and rat TIAM1, a Rho GEF regulating cell adhesion and migration via RAC1/CDC42 signaling. Ideal for cancer research.

[Learn More](#)

Parameter	Value
Product Name	TIAM1 Rabbit pAb
Remarks/Alias	NEDLDS; TIAM-1; TIAM1
Species	Human
GeneID (Human)	7074
GeneID	7074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1392-1591 of human TIAM1 (NP_003244.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13009
Protein Weight	178kDa
Shipping	Ice bag