


Tim-3/Havcr2 Rabbit Pab For Wb, Elisa - Q8Tdq0

Item Number: CM0016220



Introduction

TIM-3/HAVCR2 rabbit polyclonal antibody validated for WB and ELISA, reactive across human, mouse, and rat. Ideal for studying immune checkpoint function in T-cell regulation and cancer immunotherapy research.

[Learn More](#)

Parameter	Value
Product Name	TIM-3/HAVCR2 Rabbit pAb
Remarks/Alias	TIM3; CD366; KIM-3; SPTCL; TIMD3; Tim-3; TIMD-3; HAVcr-2; TIM-3/HAVCR2
Species	Human
GeneID (Human)	84868
GeneID	84868
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-202 of human TIM-3/HAVCR2 (NP_116171.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TDQ0
Protein Weight	33kDa
Shipping	Ice bag