

Timm13 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y5L4

Item Number: CM0024690

Introduction

TIMM13 rabbit polyclonal antibody validated for WB and ELISA in human, mouse, and rat. Targets mitochondrial import inner membrane translocase subunit Tim13 (11 kDa, UniProt Q9Y5L4), a chaperone essential for mitochondrial protein import.

[Learn More](#)

Parameter	Value
Product Name	TIMM13 Rabbit pAb
Remarks/Alias	ppv1; TIM13; TIM13B; TIMM13A; TIMM13B; TIMM13
Species	Human
Gene ID	26517
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human TIMM13 (NP_036590.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5L4
Protein Weight	11kDa
Shipping	Ice bag