

Tlr2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O60603

Item Number: CM0015172



Introduction

TLR2 Rabbit Polyclonal Antibody (pAb) validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat TLR2. Ideal for innate immunity and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	TLR2 Rabbit pAb
Remarks/Alias	TIL4; CD282; TLR2
Species	Human
GeneID (Human)	7097
GeneID	7097
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of human TLR2 (NP_003255.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O60603
Protein Weight	90kDa
Shipping	Ice bag