

Tlr4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 000206

Item Number: CM0015266



Introduction

TLR4 Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with Human and Mouse. Targets Toll-like receptor 4 involved in innate immune responses to LPS and DAMPs. Suitable for inflammation and infection research.

[Learn More](#)

Parameter	Description
Product Name	TLR4 Rabbit pAb
Remarks/Alias	TOLL; CD284; TLR-4; ARMD10; TLR4
Species	Human
Gene ID (Human)	7099
Gene ID	7099
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-228 of human TLR4 (NP_612564.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	O00206
Protein Weight	96kDa
Shipping	Ice bag