

Tm4Sf19 Rabbit Pab - Q96Dz7

Item Number: CM0021415

Introduction

TM4SF19 Rabbit Polyclonal Antibody (OCTM4) for WB, IF-P, ELISA. Recognizes human and mouse TM4SF19, a 22 kDa lysosome membrane protein regulating V-ATPase and osteoclast maturation. Suitable for immunology and bone research.

[Learn More](#)

Parameter	Value
Product Name	TM4SF19 Rabbit pAb
Remarks/Alias	OCTM4; TM4SF19
Species	Human
GeneID (Human)	116211
GeneID	116211
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human TM4SF19 (NP_001191827.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse
SWISS	Q96DZ7
Protein Weight	22kDa
Shipping	Ice bag