

Tm9Sf1 Rabbit Polyclonal Antibody For Wb And Elisa - O15321

Item Number: CM0025164

Introduction

TM9SF1 Rabbit Polyclonal Antibody for WB and ELISA, targeting human and mouse TM9SF1, an autophagy-related transmembrane protein. Ideal for research use in autophagy and membrane protein studies.

[Learn More](#)

Parameter	Value
Product Name	TM9SF1 Rabbit pAb
Remarks/Alias	MP70; HMP70; TM9SF1
Species	Human
GeneID	10548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-230 of human TM9SF1 (NP_006396.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15321
Protein Weight	69kDa
Shipping	Ice bag