

Tmem165 Rabbit Polyclonal Antibody For Wb And Elisa - Q9Hc07

Item Number: CM0022167

Introduction

Rabbit polyclonal antibody against human TMEM165, a Golgi divalent cation antiporter. Validated for WB and ELISA; ideal for glycosylation disorder and ion transport research.

[Learn More](#)

Parameter	Specification
Product Name	TMEM165 Rabbit pAb
Remarks/Alias	FT27; GDT1; CDG2K; TPARG; TMPT27; SLC64A1; TMEM165
Species	Human
GeneID (Human)	55858
GeneID	55858
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 172-230 of human TMEM165 (NP_060945.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9HC07
Protein Weight	35kDa
Shipping	Ice bag