

Tmem184A Rabbit Polyclonal Antibody For Wb And Elisa - Q6Zmb5

Item Number: CM0030846

Introduction

TMEM184A Rabbit Polyclonal Antibody, suitable for Western blot and ELISA, detects human TMEM184A (SDMG1) and cross-reacts with mouse and rat. Heparin receptor involved in vesicle transport. Ideal for vascular biology studies.

[Learn More](#)

Parameter	Value
Product Name	TMEM184A Rabbit pAb
Remarks/Alias	SDMG1; TMEM184A
Species	Human
Gene ID (Human)	202915
Gene ID	202915
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 324-413 of human TMEM184A (NP_001091089.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q6ZMB5
Protein Weight	46kDa
Shipping	Ice bag