

Tmem43 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Btv4

Item Number: CM0025534

Introduction

Rabbit polyclonal antibody against human TMEM43 (LUMA). Validated for Western blot and ELISA; cross-reacts with mouse. Suitable for studies on nuclear envelope structure, innate immune signaling via cGAS-STING, and NF-kappa-B activation. Detect a 45kDa band.

[Learn More](#)

Parameter	Value
Product Name	TMEM43 Rabbit pAb
Remarks/Alias	LUMA; ARVC5; ARVD5; AUNA3; EDMD7; AUNA2; TMEM43
Species	Human
GeneID (Human)	79188
GeneID	79188
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-310 of human TMEM43 (NP_077310.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BTV4
Protein Weight	45kDa
Shipping	Ice bag