

## Tmem44 Rabbit Pab - Q2T9K0

Item Number: CM0019917

### Introduction

TMEM44 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA, with cross-reactivity to mouse and rat. Suitable for research on human transmembrane protein 44, a membrane phosphoprotein with alternative splicing isoforms.

[Learn More](#)

| Parameter        | Value                                                                                                                  |
|------------------|------------------------------------------------------------------------------------------------------------------------|
| Product Name     | TMEM44 Rabbit pAb                                                                                                      |
| Remarks/Alias    | TMEM44                                                                                                                 |
| Species          | Human                                                                                                                  |
| GeneID (Human)   | 93109                                                                                                                  |
| GeneID           | 93109                                                                                                                  |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 110-260 of human TMEM44 (NP_001159777.1) |
| Source           | Rabbit                                                                                                                 |
| Category         | Polyclonal Antibodies                                                                                                  |
| Application      | WB, IF/ICC, ELISA                                                                                                      |
| Cross Reactivity | Mouse, Rat                                                                                                             |
| SWISS            | Q2T9K0                                                                                                                 |
| Protein Weight   | 52kDa                                                                                                                  |
| Shipping         | Ice bag                                                                                                                |