

Tmem55A Rabbit Pab For Wb, Elisa - Q8N4L2

Item Number: CM0020409

Introduction

Rabbit polyclonal antibody against human TMEM55A (PIP4P2), a type 2 phosphatidylinositol 4,5-bisphosphate 4-phosphatase involved in lipid metabolism. Validated for WB and ELISA, cross-reacts with mouse.

[Learn More](#)

Parameter	Value
Product Name	TMEM55A Rabbit pAb
Remarks/Alias	TMEM55A
Species	Human
GeneID (Human)	55529
GeneID	55529
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 130-195 of human TMEM55A (NP_061180.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8N4L2
Protein Weight	28kDa
Shipping	Ice bag