

Tnnt3 Rabbit Pab - P45378

Item Number: CM0015409

Introduction

Rabbit polyclonal antibody targeting human TNNT3 (fast skeletal muscle troponin T) for WB and ELISA applications. Cross-reacts with mouse and rat. Suitable for muscle contraction and disease research.

[Learn More](#)

Parameter	Value
Product Name	TNNT3 Rabbit pAb
Remarks/Alias	TNTF; DA2B2; beta-TnTF; TNNT3
Species	Human
Gene ID (Human)	7140
Gene ID	7140
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 147-256 of human TNNT3 (NP_001036246.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P45378
Protein Weight	32kDa
Shipping	Ice bag