

Tomm34 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q15785

Item Number: CM0028490

Introduction

TOMM34 Rabbit pAb validated for WB, IF/ICC, ELISA. Targets mitochondrial import receptor TOM34; reacts with human, mouse, rat. Ideal for mitochondrial protein import and chaperone function studies.

[Learn More](#)

Parameter	Value
Product Name	TOMM34 Rabbit pAb
Remarks/Alias	TOM34; URCC3; HTOM34P; TOMM34
Species	Human
GeneID (Human)	10953
GeneID	10953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human TOMM34 (NP_006800.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15785
Protein Weight	35kDa
Shipping	Ice bag